MPR Beckman-Weitlaner Retractor

MPR-2557



Description

A self-retaining retractor for spreading soft tissue and the edges of incisions to access a surgical site. The arms are hinged to allow use in anatomically challenging areas and the pronged tips ensure a secure grip on tissue.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	6.38” (27.5cm)
Arm Length (Shorter)	2.5” (6.35cm)
Arm Length (Longer)	2.63” (6.67cm)
Tip Length (Shorter)	0.5” (1.27cm)
Tip Length (Longer)	0.75” (1.91cm)
Prong Length	0.55” (1.4cm)
Prong Width	0.1” (0.25cm)
Tip Spread	3.88” (9.86cm)
Range of Rotation*	90°

*both directions

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com