
MPR Precision Stainless Steel Star Drivers



MPR7004



MPR7008



MPR7009



MPR7010



MPR7015



MPR7020



MPR7025

Description

Drivers for screws with star (Torx) fasteners, designed for optimal performance with minimal risk of cam-out. Made from stainless steel for durability, longevity, and easy cleaning/sterilization. MPR Orthopedics offers sizes ranging from T4 (0.05”/1.28mm) to T25 (0.173”/4.43mm), with more available upon request (*see FOR FURTHER INFORMATION*).

Made in USA.

Product Specifications

#	LENGTH	TIP SIZE
MPR7004	4.5” (11.43cm)	T4
MPR7008	5” (12.7cm)	T8
MPR7009	5” (12.7cm)	T9
MPR7010	6” (15.24cm)	T10

#	LENGTH	TIP SIZE
MPR7015	7” (17.78cm)	T15
MPR7020	7” (17.78cm)	T20
MPR7025	8” (20.32cm)	T25

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com