
MPR Capsule Retractors



MPR1054-01 MPR1054-02

Description

Retractors used during hip orthopedic surgeries to retract the acetabular capsule and expose the acetabulum. Features include rounded tips for easy tissue retraction while protecting the patient, a 130°-bent design for improved access, and a hole in the handle for attaching weights. MPR Orthopedics offers two different capsule retractors, for each hip.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Product Specifications

#	DESC. LENGTH		BLADE LENGTH	BLADE WIDTH	PRONG WIDTH	TIP WIDTH
MPR1054-01	Left	12” (30.5cm)	1.38” (3.5cm)	0.79” (2cm)	0.35” (0.9cm)	0.18” (0.45cm)
MPR1054-02	Right	12” (30.5cm)	1.38” (3.5cm)	0.79” (2cm)	0.35” (0.9cm)	0.18” (0.45cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com