MPR Narrow Screw Removal Pliers

MPR398.651



Description

Pliers to assist in removing the shafts of broken screws. The narrow jaws enhance precision and enable use in tight spaces, and feature ridged points for a better grasp on threaded shafts, and the ratchet-lock mechanism maintains a tight grip to make screw removal easier. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	7.78″ (19.75cm)
Jaw Length	0.63” (1.59cm)
Jaw Width	0.2” (0.5cm)
Point Width	0.09” (0.22cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com