
MPR Femoral/Tibial Knee Claw

MPR-00-510-486-01



Description

An extractor for removing the tibial and femoral components of a knee implant. Its claw-shaped design firmly grasps under the implant, and is secured using the screw on the end, before a slap hammer is attached to apply the extraction force. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Compatible with the Zimmer Biomet® Exodus™ Knee Revision System.

Made in USA.

Product Specifications

Minimum Length	3.38″ (8.57cm)
Maximum Length	5.5″ (13.97cm)
Minimum Point Spread	0.75″ (1.91cm)
Maximum Point Spread	2.88” (7.3cm)
Screw Size (Torx)	T40
Screw Size (Hex)	4.5mm

Product Compatibility

#	DESCRIPTION
MPR831	MPR Slap Hammer
MPR939	MPR Small Slap Hammer

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

Zimmer Biomet and Exodus are trademarks of Zimmer Biomet Holdings, Inc.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com