MPR Neck Resect Retractor

MPR-31-555516



Description

A retractor for use during Minimally Invasive Hip procedures to pull back the soft tissue surrounding the femoral neck and expose it for resection. Its blunt tip helps protect the tissue from damage while maintaining a high level of precision, and the curved blade follows the natural contours of the hip's anatomy. Features a flat handle for an optimal degree of leverage per amount of applied force, and a hole for attaching weights.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	9.87" (25.07cm)
Handle Length	3" (7.62cm)
Handle Width	1.6" (4.06cm)
Blade Length	2.3" (5.84cm)
Tip Width	0.23" (0.57cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com