MPR Toboggan Limb Guards



MPR177931



MPR177932



MPR177933

Description

Curved clear acrylic guards for preventing a patient's arm or leg from slipping off of the OR table, and protecting the patient from being leaned on or pinched during procedures. Their durable acrylic construction is designed to stand up to long-term use and is fully radiolucent, while the curved form factor allows the surgeon to lean over the patient with minimal obstrruction, and protects IV lines from accidental interference. Designed to slide right under the OR table pad, held in place by the weight of the patient.

Can be used with or without the corresponding gel pad (*see ACCESSORIES*). Made in USA.

Product Specifications

#	SIZE	HEIGHT	WIDTH	LENGTH
MPR177931	Small	6″ (15.24cm)	10″ (25.4cm)	10.25″ (26.04cm)
MPR177932	Medium	8″ (20.32cm)	16″ (40.64cm)	10″ (25.4cm)
MPR177933	Large	10.5″ (26.67cm)	15.75″ (40.01cm)	10.25″ (26.04cm)

Accessories

#	DESCRIPTION	COMPATIBILITY
MPR177930-PAD	MPR Toboggan Limb Guard Gel Pad	MPR177931, MPR177932, MPR177933

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com