MPR Posterior Glenoid Neck Retractor

MPR35000-10



Description

A retractor used during shoulder procedures utilizing a posterior approach to retract soft tissue around the glenohumeral joint and expose the glenoid neck. Its hooked handle is specially designed to allow for one-finger retraction, and the curved blade and toothed tip enable access behind the glenoid, and firmly engage with the tissue for easy manipulation.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Made in USA.

Product Specifications

Overall Length	9” (23.5cm)
Handle Length	2.5” (6.35cm)
Blade Length	2” (5.08cm)
Blade/Tip Width	1.25” (3.18cm)
Tip Length	0.75” (1.91cm)
Bend	50°

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com