MPR Vein Retractors



MPR4021062900



MPR4021082900

Description

Retractors used during a variety of procedures to retract and protect veins and other soft tissue. Feature long 90° angled blades for deep access, low-profile tapered shafts for reducing clutter and obstruction around the operating field, and fenestrated exclamation-point-shaped handles ensure a secure grip and allow for attaching weights.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Product Specifications

#	LENGTH	HANDLE LENGTH	BLADE LENGTH	BLADE WIDTH
MPR4021062900	7.96″ (20.23cm)	4.45″ (11.3cm)	5.54″ (14.08cm)	0.79″ (2cm)
MPR4021082900	7.96″ (20.23cm)	4.45″ (11.3cm)	8.13″ (20.66cm)	0.79″ (2cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com