MPR Golf Club Retractor

MPR405891



Description

A retractor used for pulling back the edges of incisions and soft tissue retraction during a variety of orthopedic and other surgical procedures. The flat, 90° bent blade provides deep access and broad retraction capabilities and features a hooked tip for a firm hold on soft tissue while maintaining a high degree of clearance, keeping the surgical field free of clutter.

Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Product Specifications

Overall Length	9” (22.86cm)
Handle Length	3.75” (9.53cm)
Blade Length	2.17” (5.5cm)
Blade Width	0.87” (2.2cm)

For Further Information

If further information on this product is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com