MPR Posterior Acetabular Retractors



MPR937500044



MPR937500045

Description

Retractors for use in hip procedures using a posterior approach to assist in acetabular exposure. Their sharp 90° angled design and flat blades help to ensure a high degree of leverage and the thin construction and small blades easily slip between anatomical structures. Feature two holes in the handles for attaching weights. MPR Orthopedics offers two Posterior Acetabular Retractors, with different working end lengths to suit your preferences or the needs of the procedure.

All MPR Orthopedics retractors are designed for optimal performance with minimal risk of trauma. Made from solid stainless steel for durability, longevity, and easy cleaning/sterilization.

Product Specifications

#	LENGTH	WORKING END	BLADE L.	BLADE W.	TIP L.
MPR937500044	10.2″ (25.9cm)	5.47″ (13.9cm)	1.38″ (3.5cm)	0.8″ (2.03cm)	0.53″ (1.35cm)
MPR937500045	10.2″ (25.9cm)	7.48″ (19cm)	1.38″ (3.5cm)	0.8″ (2.03cm)	0.53″ (1.35cm)

For Further Information

If further information on these products is required, please contact MPR Orthopedics Customer Service at 800.524.5194 in the United States, or your authorized representative, or visit the MPR Orthopedics website at www.mprortho.com.

MPR Orthopedics • 917 Lone Oak Rd, Suite 1000 • Eagan, MN 55121-2266
Phone: 800.524.5194 • Fax: 952.882.1389 • www.mprortho.com